BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

IMPORTANT: These instructions for use are designed to assist in using the Theken Spine Vu L-PoD and Vu ePoD intervertebral body fusion devices. It is not a reference for surgical techniques.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician. **Single Use!**

As with all orthopedic implants, implants and implanted system components should **never** be reused under any circumstances. An implant once used should not be reused and should be discarded. Even though it may appear undamaged, it may have small defects or internal stress patterns, which may lead to failure.

DESCRIPTION: The Vu L-PoD and Vu ePoD intervertebral body fusion devices are comprised of polymer (PEEK OPTIMA® LT) concave cages. The PEEK cages include pins (or radiopaque markers) for radiological evaluation. All pins are press fit into the PEEK material. The pins are placed in the implant on each end of the PEEK cages to allow easier radiological verification of the position and orientation of the radiocut PEEK cage.

The cages include larger central, vertical graft windows, which may be packed with autogenous bone graft material prior to implantation.

The Vu L-PoD and Vu ePoD components are manufactured from Polyetheretherketone (PEEK-OPTIMA® LT) polymer obtained from Invivo, Inc. per ASTM F-2026. Tantalum ASTM F-460 is used for the radiopaque markers.

Theken Spine warrants that these devices are fabricated from the foregoing material specification. No other warranties, expressed or implied, are made.

INDICATIONS:

When used as an intervertebral body fusion device, the Vu ePoD and Vu LPoD Intervertebral Body Fusion Devices are indicated for use as an adjunct to fusion in patients with degenerative spondylolisthesis (DO) at one or two levels from L2 to S1. These DOO patients may have up to Grade 1 spondylolisthesis at the involved level(s). The device is indicated for use with autograft only. The Vu ePoD and Vu LPoD Intervertebral Body Fusion Devices are intended for use with supplemental fixation such as the Coral Spinal System or the BodyForm Degenerative disc disease (DDD) is defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment.

CONTRAINDICATIONS:

Contraindications may be relative or absolute. The choice of a particular device must be carefully weighed against the patient's overall evaluation. Circumstances listed below may reduce the chance of a successful outcome.

- Acute or chronic infectious diseases of any etiology and localization
- Morbidity obesity
- Signs of local inflammation
- Fever or leukocytosis
- Metal/polymer sensitivity/allergies to the implant materials
- Medial or surgical conditions, which would preclude the potential, benefit of spinal implant surgery, such as the elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or a marked left shift in the WBC differential count
- Grossly distorted anatomy due to congenital abnormalities
- Rapid joint disease, bone absorption, osteopenia, and/or osteoporosis (osteoporosis is a relative contraindication since the condition may limit the degree of obtainable correction, the amount of mechanical fixation, and/or the quality of the bone graft)
- Any case not needing a bone graft and fusion or where fracture healing is not required
- Any case requiring the mixing of metals from different components
- Patients having inadequate tissue coverage over the operative site or where there is inadequate bone stock, bone quality, or anatomical definition
- Unsuitable or insufficient bone support, bone immaturity
- Any case not described in the indication
- A patient unwilling to cooperate with the postoperative instructions
- The patient's activity level, mental condition, or occupation.
- Any time implant utilization would interfere with anatomical structures or expected physiological performance
- Prior fusion at the level(s) to be treated.

These contraindications can be relative or absolute and must be taken into account by the physician when making his decision. The above list is not exhaustive.

POSSIBLE ADVERSE EVENTS:

While the expected life of spinal implant components is difficult to estimate, it is finite. These components are made of foreign materials that are placed within the body for the potential to cause injury and reduction in quality of life. However, because many biological, mechanical and physicochemical factors that affect these devices cannot be evaluated *in vivo*, the components cannot be expected to indefinitely withstand the activity level and loads of normal healthy bone. Patient should be advised of possible adverse events. This includes, but is not limited to:

- Bending, disassembly or fracture of any or all implant components
- Failure fracture of spinal fixation devices.
- Pain, discomfort, or abnormal sensations due to the presence of the device.
- Pressure on skin from components where inadequate tissue coverage exists over the implant, with the potential extension through the skin.
- Local weak requiring surgical repair. This risk is related to the surgical procedure. The intended use of the device does not require it to be close to the dura.
- Cessation of growth of the fused portion of the spine.
- Loss of proper spinal curvature, correction, height and/or reduction.
- Delayed Union or Nonunion: Internal fixation appliances are load-sharing devices, which are used to obtain alignment until normal healing occurs. In the event that healing is delayed, does not occur, or failure to immobilize occurs, the implant will be subject to excessive and repeated stresses, which can eventually cause loosening, bending or fatigue fracture.
- Loosening of spinal fixation implants can occur. Early mechanical loosening may result from inadequate initial fixation, latent infection, premature loading of the prosthesis or trauma. Late loosening may result from trauma, infection, biological complications or mechanical problems, with the subsequent possibility of bone erosion, migration and/or fracture.
- Surgical spine procedures may be associated with vascular and neural complications such as arterial injury or mechanical compromise, cord contusion and damage, peripheral nerve compromise and damage, including but not limited to peripheral paralysis, sensory deficits, vascular disorders, loss or disturbance of bladder and bowel functions.
- Serious complications may be associated with any surgery. These complications include, but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; burnitis, hemorrhage, myocardial infarction, paralysis or death.
- Neurological, vascular, or soft tissue damage due directly to the unstable nature of the fracture, or to surgical trauma.
- May increase biomechanical stress on adjacent levels.
- Inappropriate or improper surgical placement of this device may cause distraction or stress shielding of the graft or fusion mass. This may contribute to failure of an adequate fusion mass to form.
- Decrease in bone density due to stress shielding.
- Intraoperative fissure, fracture, or perforation of the spine can occur due to implantation of the components.
- Postoperative fracture of autogenous bone graft, the intervertebral body, pedicle, and/or lateral mass above and/or below the level of surgery can occur due to trauma, the presence of defects, or poor bone stock.

Adverse effects may necessitate re-operation or revision.

WARNINGS AND PRECAUTIONS:

Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without a previous surgery.

The safety and effectiveness of spinal systems have been established only for spinal conditions with significant mechanical instability or deformity requiring fusion with instrumentation. These conditions are significant mechanical instability or deformity of the thoracic and lumbar spine secondary to degenerative spondylolisthesis with objective evidence of neurological dysfunction, discitis, osteomyelitis, scoliosis, kyphosis, spinal tumor and failed previous fusion (pseudarthrosis). The safety and effectiveness of these devices for any other condition is unknown.

Only experienced spinal surgeons should perform the implantation of spinal systems with specific training in the use of this spinal system because this is a technically demanding procedure presenting a risk of serious injury to the patient. Based on the fatigue testing results, the physician/surgeon should consider the impact of implantation, patient weight, patient activity level, other patient conditions, etc, which may lead to failure of the system.

CLEANING AND DECONTAMINATION:

All instruments and implants that have been previously taken into a sterile surgical field must be first decontaminated and cleaned using established hospital methods before sterilization and reintroduction into the sterile surgical field.

- Remove all gross visible soil with a damp gauze pad or wipe.
- Prepare an enzymatic cleaning solution per manufacturer's instructions. Immerse the instruments in the cleaning solution and activate the graspers so the enzymatic cleaner contacts all mated surfaces.
- Sonicate the instruments while immersed in the cleaning solution for 15 minutes.
- Transfer the instruments to fresh enzymatic cleaning solution. Thoroughly scrub all instruments with a soft bristle cleaning brush while immersed in the enzymatic cleaning solution. Be sure that thorough scrubbing also includes any lumens with an appropriate size brush.
- Thoroughly rinse all instruments with warm running water and dry with a clean cloth and/or allow to air dry.

Note: Certain cleaning solutions such as those containing caustic soda, formalin, glutaraldehyde, bleach and/or other alkaline cleaners may damage some devices, particularly instruments; these solutions should not be used.

STERILIZATION:Unless marked sterile and clearly labeled as such, the Vu L-PoD and Vu ePoD components are provided clean but not sterile and must be sterilized prior to use. The individual products are recommended to be steam sterilized by the hospital in a gravity-displacement autoclave cycle at 270°F (132°C) for a minimum of 3 minutes exposure time with a 1 minute drying time.

The fully loaded implant and instrument tray is recommended to be steam sterilized by the hospital using an FDA cleared wrap in a gravity-displacement autoclave cycle at 270°F (132°C) for a minimum of 15 minutes exposure time with a 30 minute drying time. These sterilization recommendations follow the guidelines for sterilization per ASTM A779. Remove all packaging materials prior to sterilization. Use only sterile products in the operative field.

IMPLANT SELECTION:

The choice of proper shape, size and design of the implant for each patient is crucial to the success of the surgery. The surgeon is responsible for this choice, which depends on each patient. Once implanted, the implants are subjected to stresses and strains. These repeated stresses on the implants should be taken into consideration by the surgeon at the time of the choice of the implant, during implantation as well as in the postoperative follow-up period. Indeed, the stresses and strains on the implants may cause fatigue or fracture or deformation of the implants, before the autogenous bone graft has become completely consolidated. This may result in further side effects or necessitate the early removal of the osteosynthesis device. Improper selection, placement, positioning and fixation of these devices may result in unusual stress conditions reducing the service life of the implant.

The surgeon is to be thoroughly familiar with the surgical procedure, instruments and implant components prior to performing surgery. Refer to the Theken Spine surgical protocols for additional procedural information. Periodic follow-up is recommended to monitor the position and state of the implants, as well as the condition of the adjoining bone.

PREOPERATIVE:

Only patients that meet the criteria described in the indication should be selected. Patient conditions and/or pre-dispositions such as those addressed in the aforementioned contraindications should be avoided.

The surgeon must ensure that all necessary implants and instruments are on hand prior to surgery.

Care should be used in the handling and storage of the implant components. The implants should not be scratched or otherwise damaged. Implants and instruments should be protected during storage especially from corrosive environments.

They should be carefully unpacked and inspected for damage prior to use.

The type of construct to be assembled for the case should be determined prior to beginning surgery. An adequate surgery. An adequate surgery. An adequate surgery. An adequate surgery. Including sizes larger and smaller than those that are expected for use.

The correct handling of the implant components is extremely important. Contouring of the devices is to be avoided.

Because mechanical parts are involved, the surgeon should be familiar with the various components before using the equipment and should personally assemble the devices to verify that all parts and necessary instruments are present before the surgery begins.

INTRAOPERATIVE:

At all times, extreme caution should be used around the spinal cord and nerve roots. Damage to nerves may occur resulting in loss of neurological functions.

The implant surfaces should not be scratched or notched since such actions may reduce the functional strength of the construct.

Autogenous bone grafts must be placed in the area to be fused and the graft must be in contact with visible bone.

POSTOPERATIVE:

The physician's postoperative directions and warnings to the patient and the corresponding patient compliance are extremely important.

Detailed instructions on the use and limitations of the device should be given to the patient. If partial weight-bearing is recommended or required prior to firm bony union, the patient must be warned that bending, loosening or breakage of the components are complications which can occur as a result of excessive or early weight-bearing or excessive muscular activity. Bending, loosening, or breakage of an internal fixation device during postoperative rehabilitation may be increased if the patient is active, or if the patient is debilitated, demented, or otherwise unable to use bracing or other such weight supporting devices. The patient should be warned to avoid falls or sudden jolts in spinal position. To allow maximum chances for a successful surgical result, the patient or device should not be exposed to mechanical vibrations that may loosen the device construct. The patient should be warned of this possibility and instructed to limit and restrict physical activities, especially lifting, twisting motions and any type of sport participation. The patient should be advised not to smoke or consume alcohol during the autogenous bone graft healing process.

The patient should be advised of their inability to bend at the point of spinal fusion and taught to compensate for this permanent physical restriction in body motion.

Dural leak requiring surgical repair. This risk is related to the surgical procedure. The intended use of the device does not require it to be close to the dura. If a nonunion develops or if the components loosen, bend, and/or break, the device(s) should be revised and/or removed immediately before serious injury occurs. Failure to immobilize a delayed nonunion of bone will result in excessive and repeated stresses on the device(s). By the mechanism of fatigue these stresses can cause eventual bending, loosening, or breakage of the device(s).

Expanded surgical implants must never be reused. Adequately instruct the patient in the appropriate postoperative care. The patient's ability and willingness to follow instructions is one of the most important aspects of successful healing.

REMOVAL OF IMPLANT AFTER HEALING:

Particularly in young active patients, implants may loosen, fracture, migrate, increase the risk of infection, cause pain, or stress shield bone – even after normal healing. The surgeon should consider the risks and benefits when deciding whether or not to remove an implant. Implant removal should be followed by careful postoperative management to avoid re-fractures. If the patient is older and risks are low a activity level, the surgeon may elect not to remove the implant in order to eliminate the risks of an anterior surgery.

PACKAGING:

The implants are delivered in packages; these must be intact at the time of receipt. The systems are sometimes supplied as a complete set; implants and instruments are arranged on trays and placed in specially designed storage boxes.

PRODUCT COMPLAINTS:

Any health professional having a complaint or grounds for dissatisfaction relating to the quality of the product, its identity, its durability, its reliability, its safety, its effectiveness and/or its performance, should notify Theken Spine or its representative. Moreover, if a device has malfunctioned, Theken Spine or its representative must be advised immediately.

If a Theken Spine product has ever worked improperly and could have caused or contributed to the serious injury or death of a patient, the distributor or Theken Spine Representative must be informed as soon as possible by telephone, fax or in writing.

For all complaints, please include the device name and reference along with the lot number of the component(s), your name and address and an exhaustive description of the event to help Theken Spine understand the causes of the complaint.

PRODUCT INFORMATION DISCLOSURE:

THEKEN SPINE HAS EXERCISED REASONABLE CARE IN THE SELECTION OF MATERIALS AND THE MANUFACTURE OF THESE PRODUCTS. THEKEN SPINE DISCLAIMS ALL WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THEKEN SPINE SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY, ARISING FROM THE USE OF THIS PRODUCT. NEITHER THEKEN SPINE ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS. THEKEN SPINE INTENDS THAT THIS DEVICE SHOULD BE USED ONLY BY PHYSICIANS HAVING RECEIVED PROPER TRAINING IN THE USE OF THE DEVICE.

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